Your Partner
We can do it

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


FORMATION OF A THYMINE-THYMINE DIMER

FORMATION OF A THYMINE-CYTOSINE DIMER

FORMATION OF A 6-4 LESION

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펄스UV살균 - SUMMARY

■ EFFECTIVE INACTIVATION

- Achieve USP 10-6 SAL
- No microorganism found resistant to Pulsed UV
- No repair of nucleic acids

■ SELECTIVE BIOMOLECULAR DESTRUCTION

- Destruction of DNA/RNA
- High antigen recovery
- No significant heat effect on product

■ RAPID & CONTROLLED PROCESS

- Single step process
- Treatment in minutes (per liter)
- Broad spectrum fluence monitoring

■ SAFE


- No toxic chemicals or ionizing irradiation

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Sterilization of Static Liquids with Pulsed UV


Organism	Inoculation (in 20 mL)	Sterility Test Result ¹		
		Glucose	Saline	WFI
<i>Bacillus subtilis</i> var. niger spores	4.8 x 10	Pass	Pass	Pass
<i>Bacillus pumilus</i> spores	3.0 x 10	Pass	Pass	Pass
<i>Bacillus stearothermophilus</i> spores	4.4 x 10	Pass	Pass	Pass
<i>Clostridium sporogenes</i> spores	2.4 x 10	Pass	Pass	Pass
<i>Aspergillus niger</i> spores	1.4 x 10	Pass	Pass	Pass
<i>Candida albicans</i>	3.0 x 10	Pass	Pass	Pass
<i>Deinococcus radiodurans</i>	2.2 x 10	Pass	Pass	Pass
<i>Staphylococcus aureus</i>	5.6 x 10	Pass	Pass	Pass
<i>Enterococcus faecalis</i>	2.0 x 10	Pass	Pass	Pass
<i>Escherichia coli</i>	3.4 x 10	Pass	Pass	Pass
<i>Salmonella choleraesuis</i>	1.8 x 10	Pass	Pass	Pass
<i>Pseudomonas aeruginosa</i>	2.8 x 10	Pass	Pass	Pass

¹ Treatment level was approximately 4.5 J/cm² per flash, with 2 flashes. 20 replications for each liquid/microbe combination. No dark or light enzymatic repair of nucleic acids after 28 days of incubation.

Reference: Furukawa et al., Japanese PDA, 1999

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Pulse UV Demonstrates Complete Inactivation of a Cancer Causing virus with High Antigen Recovery

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


펄스UV살균 시스템 REFERENCE LIST

COMPANY NAME	SYSTEMS SOLD	APPLICATION
Advance Therapeutic Technologies	RS-3000; RC-742	Medical Device sanitization
Ainia Centro Tecnológico (Spain)	RS-3000	Packaging material (PET, multi-layer paper; aluminum) for Yogurts
Alabama A&M University	RS-3000	Food Investigation
Alcon	RC-742, RS-3000	Decontamination: Poly Bags
Biogentis, Inc (Canada)	RS-3000M-4	Bio Fill Package System online: to sterilize a medical device: ampoule.
Branchy Technology, Co (Taiwan)	RC-747	Unknown
California Poly State University	RC-847	Decontamination: Food
Cooper Vision (UK)	RS-3000; RC-747	Decontamination: Contact Lens
Cornell University	RS-3000	Decontamination: Food
Cryovac	RC-747	Food Investigation
Dayla Liquid Packaging, Limited (UK)	RS-3000	Packaging Machine : Safe Fill
EADS	RC-747	Unknown
Edelweiss	RC-847	Shelf life of pies without preservatives

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펄스UV살균 시스템 REFERENCE LIST

COMPANY NAME	SYSTEMS SOLD	APPLICATION
Entemann's	RS-3000	Increase shelf life of bakery foods
Facultad De Cienias	RC-747	Food
Fogg Filler	RC-801, RS-3000	Decontamination: Bottle Caps
Fort Valley State	RS-3000	Food Investigation
Georgia Tech University	RS-3000	Food Investigation
Iowa State University	RS-3000	Food Investigation
Little Pearl	RC 747	Increase shelf life of fish eggs
McGill University	RS-3000	Food Investigation
Millenia Healthcare Corporation	RC-742	Food
Nature Sunshine Product	RC-747	Food
Novartis	RC-747	Spore Reduction
Pepsico, Inc.	RC-801	Bottle Cap
ReVision Optics	RC-742	Implanted Eye Lens sterilization
Rutgers University	RS-3000	Food Investigation

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펄스UV살균 시스템 REFERENCE LIST

COMPANY NAME	SYSTEMS SOLD	APPLICATION
State of FL	RS-3000	Food Investigation
Steris	RC-747	Medical Device sterilization
University of Salerno	RC-801	Food Investigation
University of Alaska	RC-747	Food Investigation
University of Auckland	RC-500	Food Investigation
University of Manitoba	RS-3000	Food Investigation
USDA	RS-3000	Food Investigation
Etc.		

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- AINIA Centro Tecnologico, Laboratorio de Bioensayos, Valencia, Spain
- Alabama A&M
- American Red Cross Holland Laboratories
- Collage of Food Science –FugianAgriculture and Forestry University, China
- Contamination Zero
- Cornell University
- Florida Dept of Citrus
- Fort Valley State University
- Georgia State University
- Institute of Food Technology –Taiwan
- Nissan Pharmaceuticals
- Nutramax Corp
- Pennsylvania State University
- Pure Pulse Technologies
- Rutgers University
- University of Massachusetts
- etc.

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Xenon Corporation